



ihcDirect™ Tools

Introducing the Novodiox ihcDirect™ Slide Warmer System and ihcDirect™ Slide Pan Assembly.



The ihcDirect Slide Warmer and Slide Pan Assembly support the Novodiox test kits by consolidating test workflow into a small area. The Warmer and Pan Assembly supports quality testing by helping users to obtain consistent results via precise temperature settings.

ihcDirect Slide Pan Assembly features include;

- Supports ihcDirect IHC assays.
- Compact design.
- Allows users to simultaneously manage up to 9 slides at a time.
- Keeps wash buffer separated from users hands by providing sufficient capacity to manage used wash buffer from a complete run.



Product Specifications

ihcDirect Slide Warmer	M81500-110
Overall size	10¼" x 8" x 3½" (260 x 200 x 85mm)
Surface size	10" x 7" (25.4 x 17.8 cm)
Capacity	Approx. 23 slides
Power	100 Watts
Net Weight	7 lbs.
Voltage	110V, 60Hz or 220V, 50Hz (M81500-220)
Temperature	Room temperature to 75° C; ±2° C
Quality Assurance	ISO 9002

ihcDirect Slide Pan Assembly	M81050
Overall Size	10 ⁷ / ₁₆ " x 6" x 1 ¹ / ₄ " (265 x 153 x 31mm)
Wash Fluid Capacity	Outer well – 40mL Center Well 125mL
Slide Size Dimension	25.5-26.5 +/- 0.5mm

ihcDirect Slide Warmer made by C & A Scientific

Janet Yang
 Director of Marketing & Business Development
janetyang@novodiox.com
 Mobil: 650-759-4500